

# ProveBlue™

## Typical Analysis

Item	Typical analysis of Provepharm Methylene Blue	Limits	
		Ph.Eur	USP*
Al	< 2	< 100	-
As	Not detected	-	< 6
Cd	< 0,01	< 1	-
Cr	< 3	< 100	-
Cu	< 3	< 300	< 200
Sn	< 0,2	< 10	-
Fe	< 3	< 200	-
Mn	< 0,5	< 10	-
Hg	< 0,5	< 1	-
Mo	< 1	< 10	-
Ni	< 2	< 10	-
Pb	Not detected	< 10	-
Zn	< 2	< 100	Absence